Bench aid handbook

for

Microscopic Diagnosis

of

Malaria

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Complied by Robin Wales

With grateful acknowledgement to;

Royal Perth Hospital Malaria Resource. Copyright State of Western Australia

Shoklo Malaria Research Unit

CDC Public Image Library

World Health Organisation


Tropical Health Technology and Professor Balbir Singh for P. knowlesi information.


Burmese translation courtesy Saw Marcus and U Maung Maung Win.


edited by Elizabeth Nyein Nyein New

All pictures are taken using 100 × oil immersion objective in Geimsa stained thick and thin blood films.

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Normal blood film. Red and white blood cells and a platelet. No malarial parasites.

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Microscopic examination of blood films

It is important that blood films for microscopy are well made and properly stained. Thin films must be spread evenly so that individual RBCs can easily be seen and examined in the thinner end of the film. Methanol used for thin film fixation must be pure and contain no water. The slides must be kept dry and fixed as soon as possible. Any contact with water will destroy the thin film.

Thick films can be unusable if the RBCs become fixed and are not washed away during staining. This may happen if slides are stored for too long or put in hot places. It can also happen if methanol, or the vapour from the methanol container, comes in contact with the thick film when the thin film is being fixed.

Thick films are used for diagnosis of malaria. It is only necessary to examine the thin film in order to identify the species if the thick film was positive, or to report the percentage of RBCs containing Pf parasites. This is sometimes done during treatment to check the result of medication.

Thick films should first be examined using the 10x and 40x microscope objectives. We do this to select best stained parts of the slide for examination. The best part of the thick film has about 10 White blood cells per High Power Field (HPF). It is sometimes possible to identify other things at this stage, including Microfilaria. We may also see dark brown Malaria pigment in WBCs, this is an important finding even if no parasites are found because it is a clue that malaria parasites have been in the cells.

Raise the objectives above the slide and put a small drop of immersion oil on the selected area before examining with the 100× objective.

Examine the film using a pattern as below, moving from one field to the next systematically, continue for 100 good fields. Report result as positive or neagative for malaria. If positive record species, stages and numbers of parasites.

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This diagram shows the correct way to view the thick and thin sections of a blood film.

If the thick film is unusable (fixed, washed off, etc,) but the thin film is good, examine 500 fields in the thin film. If this result is negative a fresh slide should be requested and the thick film reported from the new slide.

Most microscopes are supplied with 10× eyepieces. If possible, it is better to use 7× or 8x eyepieces as they give a clearer image with the 100× oil immersion objective. Never use eyepieces over 10×.

Reporting parasite numbers can be done in different ways. For routine reporting in thick films, with an average of 10 WBCs per HPF, use the following:

Number of Parasites	*Report as*
Less than one parasite per HPF	Scanty or Rare
1—2 parasite(s) per HPF	+ (1+)
3 –25 parasites per HPF	++ (2+)
26—60 parasites per HPF	+++ (3+)
Over 60 parasites per HPF	++++ (4+)

When thick film parasite numbers are 3+ or 4+, an accurate count is required to determine treatment of hyperparasitaemia. To do this, the thin film must be examined to express parasite density using the percentage of RBCs infected. This is done by finding a part of the thin film where RBCs are evenly distributed, and counting the number infected RBCs in fields of view containing a standard number of RBCs. For example count the number of infected RBCs in 10 fields with an average of 200 RBCs in each. This will be the number per 2,000 RBCs, divide this number by 20 to report as a percentage.

Sometimes parasite counts are made by comparing with the number of WBCs in thick films.

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
Thick film showing no malarial parasites. At least 100 fields must be examined before reporting as negative.

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
Always use new or properly cleaned slides for making blood smears. Unclean slides, poor quality Geimsa stain, or water which is too acid or alkaline can result in unsatisfactory staining.

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If the water used for diluting the Giemsa stain is slightly alkaline, the image will be too blue. If the water is slightly acid, the image will be too red. The pH of the water can be corrected using a pH 7.2 buffer made by adding buffer concentrate or buffer tablets to water used for stain dilution.

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
Developing P.vivax trophozoite. Because the water used to dilute the Giemsa stain was alkaline, the slide is very blue and Schüffner's dots are not visible..

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
P.vivax trophozoite. This is a well stained film. Note Scuffner's dots, the red cells and Pv trophozoite are pinker in this slide.

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P.ovale trophozoite. Because it is well stained the James' dots are clearly visible.

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P.ovale trophozoite. The water used to dilute the Giemsa stain was slightly acid, making the film appear very red. As a result the James' dots are only just visible.

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Even using the same stain and staining method there will be variations in the colour and quality of staining in the slides you examine. You will notice these differences in stain quality and parasite appearance in the photos in these notes. It is important to learn to recognise parasites in differently stained slides.

It is important for you to be able to recognise things which can be mistaken for a malaria parasite. These things include stain crystals, platelets and debris from dirty slides.

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There are four species of malaria primarily infecting humans. The parasites are transmitted by mosquitoes. Two of these, *Plasmodium falciparum* (Pf) and *Plasmodium vivax* (Pv) are the most common in Southeast Asia and together cause over 95% of cases of Malaria. Overall occurrence of each is about 50%, with Pf more frequent in tropical areas and Pv in more temperate places.

Pf is the most important and causes many deaths, especially in children. It also causes serious problems for pregnant women and their unborn and newborn babies. We need to identify this quickly by microscopy or RDT because special medicine must be given as soon as possible.

Pv is also an important cause of illness and loss of income.
We can diagnose malaria and also see which species is causing it by malaria microscopy. The parasites are found inside red blood cells and differences in appearance of these for each species enable us to differentiate them.

The other two species, *Plasmodium ovale* Po, and *Plasmodium malariae* Pm, have both been found in Thailand and Burma, but are not common. Like Pv they are less serious than Pf and respond well to treatment.

Plasmodium knowlesi (Pk) is a malarial parasite of monkeys which more and more commonly causing severe malarial illness in man. It is mainly found in Malaysia where large areas of forest have been cleared causing much higher

risk of mosquitoes transferring the disease to humans.

We will be using the scientific names given to the forms of malarial parasites at different stages of their growth and maturation inside red blood cells as seen on page 15.

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

Tiny single cell merozoites (pictured) on left from blood or liver schizonts enter red blood cells and become early ring form

trophozoites. Shown in the picture on the right with red nuclear dot and blue cytoplasm.

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
This diagram shows a comparison of all the Plasmodium species within their developmental stage. The thin film is shown on the top line and the thick film on the bottom line. The size and shape of parasites and infected red blood cells varies between species.

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P.falciparum


This diagram shows the developmental stages of P. falciparum. We do not usually see the stages of P. falciparum parasites from young trophozoites to mature schizonts (shown in boxes) when we examine blood films. These stages occur in capillary blood vessels of deep inner organs and in the placenta of the pregnant woman.

Diagnosis of *P. falciparum* malaria is based on finding the ring form trophozoites and gametocytes characteristic of this species in thick and thin blood films. The typical shaped gametocytes enable us to make a certain diagnosis of falciparum malaria.

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This diagram shows the appearance of the Pf parasite in the thin films on the left hand side, and thick films on the right. In this species the size of the red cells is not increased and they may have crenated (ragged) edges. Dots developing as trophozoites mature are sometimes seen in the thin films. These are called Maurer's clefts.

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Identifying features of P.falciparum


- Usually only ring trophozoites and gametocytes are seen.
- The red cell size is not enlarged.
- The red cells are often crenated (jagged edges).
- Large numbers of small ring trophozoites are often present.
- Sometimes more than one parasite is seen in the red cells. Often two nuclei are seen. They may be on the edge of the red cells.
- No regular dots like the Schüffner's dots of Pv are present in early rings, but irregular dots called Maurer's clefts may be seen in maturing trophozoites.
- The gametocytes have a typical banana shape. This is diagnostic of Pf. Malaria pigment granules may be present, mainly around nucleus.

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
Thick films of P.falciparum

Thick film showing Pf ring trophozoites and a gametocyte.

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

Thick film showing Pf growing trophozoites.

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
Thick film showing Pf ring trophozoites.

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
Thick film showing Pf ring trophozoites. Notice the difference in staining qualities of this image and the previous one.

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
Thick film showing mixed infection of Pf gametocytes and Pv trophozoites.

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
Thick film showing many Pf rings, four distinct gametocytes, and three distinct schizonts. Schizonts are small and rarely seen in the circulating blood, except in very severe malaria. They contain two to four, and occasionally up to eight merozoites.

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Thick film showing malaria pigment in two white blood cells.

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
Thick film showing a few Pf rings and many platelets.

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
Thin films of P.falciparum


Thin film showing malaria pigment in a WBC. This is an important finding and must always be reported.

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
Thin film showing platelet artefacts and a normal Pf gametocyte. It is important to be able to tell the difference between these..

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
Thin films showing a female gametocyte of Pf on the left and male gametocyte of Pf on the right.

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
Thin film showing growing Pf trophozoites with Maurer's clefts and crenated red cells.

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Thin film showing Pf rings with three on the wall of one red cell. This is called the accolé form.

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
See next page

Thin film shows heavy Pf infection with multiple rings in RBCs, double dots and crenated RBCs.

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P.vivax

This diagram shows the developmental stages of P. vivax. All stages of trophozoites, schizonts, and gametocytes can be found in thick and thin blood films.

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Identifying features of P.vivax

- Red cell size increases greatly as parasites grow and mature. The RBC is often irregular in shape due to the amoeboid nature of the parasite.
- Spreading cytoplasm becomes irregular and fragmented during growth and many Schüffner's dots are seen in well stained films.
- Few to moderate numbers of parasites are present unlike the very large numbers often seen in Pf infections.
- Occasionally more than one ring form and double nuclei are seen in

red cells, but this is not seen as commonly as it is in falciparum infections.


- The gametocytes are usually large, round or irregular, and when mature they almost completely fill the red cell. Small forms may also be present.
- Schizonts also large round or irregular, filling the RBCs when mature. They are often ruptured in thick films and reveal 16 to 24 or more merozoites.

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


Thick film showing artefacts including platelet clump. No parasites are present.

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Previous page - Thick film showing Pv gametocyte and small trophozoite.

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
Thick film showing Pv trophozoites and schizonts. Mature schizonts may contain 24 or more merozoites. Note the pigment granule in the schizont above centre.

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
Thin films of P.vivax


Thin film showing young Pv trophozoites. You can see two rings and two growing rings. The RBC is increasing in size and amoeboid shape. Note also the Schüffner's dots

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Thin film showing larger Pv trophozoite and schizont showing clear Schüffner's dots.

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
See next page


Thin film showing two Pv rings in one enlarged RBC and two amoeboid growing trophozoites with Schüffner's dots.

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
Thin film showing Pv mature schizont.

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
(1)Thin film showing Pv microgametocyte.(male)

(2)Thin film showing pv macrogametocyte.(female)

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P.ovale

This diagram shows the developmental stages of Po.This species is rarely found in Southeast Asia. All stages of trophozoites, schizonts, and gametocytes can be found in blood films. Usually no more than 2% of RBCs are infected. Occasionally more than one parasite found in RBCs.

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Identifying features of P.ovale

- As with Pv, red cell size increases as parasites grow and matures. However, the parasites are smaller than those of Pv and do not occupy the whole host RBC. Up to about 30% of infected RBCs become oval


shaped or have fimbriated (jagged) ends. As trophozoites grow, James' dots appear. They are similar to the Schüffner's dots of Pv, but more obvious. They are found in both thick and thin films.

- The gametocytes are small and usually round, with James' dots.
- The schizonts are small with up to 10 merozoites and James' dots. The RBCs are often oval with jagged ends. Only small amounts of malaria pigment may be seen, mainly in the schizonts.

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
Thick and thin films of P.ovale

Thick film showing Po schizont with 8 separated merozoites and a growing trophozoite above. The trophozoite shows a large red blood cell outline. An immature schizont is also seen with two chromatin masses. The parasites are smaller than those seen in Pv

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

Thin film showing two RBCs with growing Po trophozoites and James' dots. One RBC has two parasites.

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
Thin films showing a Po microgametocyte (male) on the left, and a Po macrogametocyte (female) on the right.

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
Thin films showing Po early trophozoite on the left, and Po growing trophozoite on the right.

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
film showing Po mature trophozoite. This can be confused with the Po gametocyte

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P. malariae

The previous diagram shows the developmental stages of P. malariae. This parasite is only sometimes found in Southeast Asia.

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Identifying features of P. malariae

- All stages of trophozoites, schizonts, and gametocytes are found in the blood films. Rarely more than 1% of RBCs are infected.
- RBCs are not enlarged and do not contain dots. Ring forms, band forms and basket forms of trophozoites may be present. Malaria pigment may be present in growing trophozoites. Mature trophozoites appear similar to gametocytes.
- The gametocytes are small, round or oval, with yellow brown malaria pigment. The nucleus usually lies to one side.
- The schizonts are small with up to 12 merozoites and malaria pigment.
 Mature schizonts often have merozoites neatly arranged around malaria pigment to give a daisy like appearance.

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

Thick films of P. malariae


Thick film showing a Pm schizont and early trophozoite.

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
Thick film showing Pm late trophozoites and developing schizont.

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previous page: Thick film showing ruptured Pm schizont with merozoites and malaria pigment

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
Thick film showing Pm gametocyte and trophozoite.

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
Thin films of P. malariae


Thin film showing Pm early ring trophozoite and schizont

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
Thin film showing Pm band form early trophozoites.

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
Thin film showing Pm late trophozoite band form.

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
Thin film showing Pm basket form trophozoite.

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
Thin film showing Pm developing schizont.

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
Thin film showing Pm mature schizont.

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P. knowlesi

This is a Malarial parasite of monkeys that can also infect humans. It is mainly in Malaysia and causes a serious and life threatening disease. In some areas it accounts for half the human malaria. Sometimes it is found as a mixed infection with Pf. This species of *Plasmodium* has been found in Thailand and could become a problem in Burma as forests are cleared and humans have closer contact with monkeys. Due to similarity in appearance it can be wrongly diagnosed as Pm In this infection many RBCs are infected and all stages may be present in smears at the same time. The blood cycle is 24hours with fever recurring daily.

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Identifying features of P. knowlesi

- RBCs are not enlarged. They may have fimbriated ends and may contain small dots called Sinton stippling.
- The trophozoites may look like Pf rings if only early forms are present. There may be more than one parasite and some with double dots in RBCs. The cytoplasm is compact and band and basket forms like those of Pm may be found.
- The schizonts resemble those of Pm with malaria pigment. They contain up to 16 merozoites.
- The gametocytes are also round with scattered malaria pigment. They tend to fill most of the Red blood cell. The overall appearance is similar to that of the Pm gametocytes.

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

Thick films of P. knowlesi

Previoius page: Thick film showing a heavy infection of Pk with trophozoites, gametocytes and schizont.

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
Thick film showing of Pk trophozoites and gametocytes.

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
Thin films of P. knowlesi


Thin film showing Pk trophozoites, some with dots and a basket form

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Thin film with Pk Gametocytes and band trophozoites

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Thin film showing Pk schizonts and a trophozoite

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